


Rabbit Anti-Mouse Cd206/Mmr Monoclonal Antibody For Fc - Q61830

Item Number: CM0011520



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD206 (macrophage mannose receptor 1). Designed for flow cytometry analysis of mouse macrophages. Specific for mouse CD206, molecular weight approximately 165 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD206/MMR mAb
Remarks/Alias	MR; CD206
Species	Mouse
GeneID (Human)	4360
GeneID	17533
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-154 of mouse CD206/MMR (NP_032651.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra),FC
Cross Reactivity	Mouse
SWISS	Q61830
Protein Weight	165 kDa
Shipping	Ice bag