


Rabbit Anti-Mouse Cd206/Mmr Monoclonal Antibody For Flow Cytometry - Q61830

Item Number: CM0009845



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD206/MMR for flow cytometry. Validated for intra- and surface staining. Ideal for macrophage and phagocytosis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD206/MMR mAb
Remarks/Alias	MR; CD206
Species	Mouse
GeneID (human)	4360
GeneID	17533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-154 of mouse CD206/MMR (NP_032651.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra), FC
Cross Reactivity	Mouse
SWISS	Q61830
Protein Weight	165 kDa
Shipping	Ice bag