


Rabbit Anti-Mouse Cd209B/Dc-Signr1 Monoclonal Antibody For Flow Cytometry - Q8Cj91

Item Number: CM0009727



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD209b/DC-SIGNR1, validated for flow cytometry. Suitable for research on dendritic cell biology, pathogen recognition, and immune receptor studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD209b/DC-SIGNR1 mAb
Remarks/Alias	OtB7; SIGNR1; mSIGNR1; DC-SIGNR1; 1810030I22Rik
Species	Mouse
GeneID (Human)	30835
GeneID	69165
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 76-325 of mouse CD209b/DC-SIGNR1 (NP_081248.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8CJ91
Protein Weight	16kDa/33kDa/36kDa/37kDa
Shipping	Ice bag