

Rabbit Anti-Mouse Cd21/Cd35 Monoclonal Antibody For Fc - P19070

Item Number: CM0011772

Introduction

Rabbit anti-Mouse CD21/CD35 monoclonal antibody validated for flow cytometry. Specific to mouse CD21/CD35 (CR1/CR2), this antibody supports B-cell complement receptor studies, including B-cell activation and IgE regulation. SWISS: P19070.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD21/CD35 mAb
Remarks/Alias	Cr1, C3DR, CD21, CD35, Cr-1, Cr-2
Species	Mouse
Gene ID (Human)	1380
Gene ID	12902
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-274 of mouse CD21/CD35 (NP_031784.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P19070
Protein Weight	113kDa
Shipping	Ice bag