

Rabbit Anti-Mouse Cd22 Monoclonal Antibody For If-P And Fc - P35329

Item Number: CM0011893

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD22 (Siglec-2). Validated for immunofluorescence (IF-P) and flow cytometry (FC) on mouse B-cells. Ideal for B-cell biology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD22 mAb
Remarks/Alias	Lyb8; Lyb-8; A530093D23
Species	Mouse
Gene ID (Human)	933
Gene ID (Mouse)	12483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-151 of mouse CD22 (NP_033975.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, FC
Cross Reactivity	Mouse
SWISS	P35329
Protein Weight	89kDa / 93kDa / 97kDa
Shipping	Ice bag