

Rabbit Anti-Mouse Cd22 Monoclonal Antibody For If-P And Fc - P35329

Item Number: CM0008662

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD22. Validated for immunofluorescence (IF-P) and flow cytometry (FC). Suitable for B-cell biology studies involving differentiation, antigen presentation, and trafficking.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD22 mAb
Remarks/Alias	Lyb8; Lyb-8; A530093D23
Species	Mouse
GeneID (Human)	933
GeneID	12483
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-151 of mouse CD22 (NP_033975.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, FC
Cross Reactivity	Mouse
SWISS	P35329
Protein Weight	89kDa/93kDa/97kDa
Shipping	Ice bag