


Rabbit Anti-Mouse Cd226/Dnam-1 Monoclonal Antibody For Fc - Q8K4F0

Item Number: CM0004997



Introduction

-conjugated rabbit monoclonal antibody specific for mouse CD226/DNAM-1 (SWISS Q8K4F0), validated for flow cytometry. CD226 is a costimulatory receptor on NK cells and CTLs, enhancing cytotoxicity against infected or tumor cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD226/DNAM-1 mAb
Remarks/Alias	Pta1; DNAM1; DNAM-1; TLISA1
Species	Mouse
Gene ID (Human)	10666
Gene ID (Mouse)	225825
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-254 of mouse CD226/DNAM-1 (NP_848802.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8K4F0
Protein Weight	38kDa
Shipping	Ice bag