


Rabbit Anti-Mouse Cd226/Dnam-1 Monoclonal Antibody () For Flow Cytometry - Q8K4F0

Item Number: CM0008532



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD226/DNAM-1. Validated for flow cytometry to study immune cell adhesion, NK and CTL cytotoxicity, and lymphocyte signaling. Reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD226/DNAM-1 mAb
Remarks/Alias	Pta1; DNAM1; DNAM-1; TLISA1
Species	Mouse
Gene ID (Human)	10666
Gene ID	225825
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-254 of mouse CD226/DNAM-1 (NP_848802.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8K4F0
Protein Weight	38kDa
Shipping	Ice bag