

## Rabbit Anti-Mouse Cd23 Monoclonal Antibody For Flow Cytometry - P20693

Item Number: CM0008573

### Introduction

-conjugated rabbit monoclonal antibody against mouse CD23 (FceRII), validated for flow cytometry. Low-affinity IgE receptor involved in B cell differentiation and IgE regulation.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD23 mAb                                                                                                            |
| Remarks/Alias    | CD23; Fce2; Fcer2; Ly-42                                                                                                              |
| Species          | Mouse                                                                                                                                 |
| Gene ID (Human)  | 2208                                                                                                                                  |
| Gene ID (Mouse)  | 14128                                                                                                                                 |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-331 of mouse CD23 (NP_038545.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | FC                                                                                                                                    |
| Cross Reactivity | Mouse                                                                                                                                 |
| SWISS            | P20693                                                                                                                                |
| Protein Weight   | 37kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |