

Rabbit Anti-Mouse Cd23 Monoclonal Antibody For Flow Cytometry - P20693

Item Number: CM0010709

Introduction

Explore this -conjugated rabbit monoclonal antibody against mouse CD23 (low-affinity IgE receptor, FcεR2). Validated for flow cytometry, ideal for B cell differentiation and IgE regulation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD23 mAb
Remarks/Alias	CD23; Fcε2; FcεR2; Ly-42
Species	Mouse
GeneID (Human)	2208
GeneID	14128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-331 of mouse CD23 (NP_038545.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P20693
Protein Weight	37KDa
Shipping	Ice bag