

Rabbit Anti-Mouse Cd23 Monoclonal Antibody For Flow Cytometry - P20693

Item Number: CM0006168

Introduction

Rabbit anti-Mouse CD23 (Fcer2) monoclonal antibody, validated for flow cytometry. Recognizes low-affinity IgE receptor involved in IgE regulation, B-cell differentiation, and macrophage activation. Ideal for allergy and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD23 mAb
Remarks/Alias	CD23; Fce2; Fcer2; Ly-42
Species	Mouse
GeneID (Human)	2208
GeneID	14128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-331 of mouse CD23 (NP_038545.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P20693
Protein Weight	37KDa
Shipping	Ice bag