

Rabbit Anti-Mouse Cd24 Mab For Fc - P24807

Item Number: CM0005224

Introduction

Rabbit anti-Mouse CD24 monoclonal antibody for flow cytometry detection of mouse CD24. Validated species reactivity: mouse. Ideal for immune cell profiling and autoimmune disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
GeneID (Human)	100133941
GeneID	12484
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag