

Rabbit Anti-Mouse Cd24 Monoclonal Antibody For Flow Cytometry - P24807

Item Number: CM0010387

Introduction

-conjugated rabbit monoclonal antibody against mouse CD24 (HSA, Ly-52, nectadrin) for flow cytometry. Recognizes mouse CD24. Suitable for immune cell differentiation and activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
Gene ID (Human)	100133941
Gene ID (Mouse)	12484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag