

Rabbit Anti-Mouse Cd24 Monoclonal Antibody For If/Fc - P24807

Item Number: CM0008668

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD24, validated for immunofluorescence and flow cytometry. Reliable tool for immune cell research involving B-cell activation and thymocyte development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD24 mAb
Alias	HSA, Cd24, Ly-52, nectadrin
Species Reactivity	Mouse
Human Gene ID	100133941
Mouse Gene ID	12484
Immunogen	Recombinant fusion protein containing amino acids 27-53 of mouse CD24 (NP_033976.1)
Source	Rabbit
Antibody Type	Monoclonal Antibodies
Applications	IF/ICC, IF-P, FC
Cross-Reactivity	Mouse
UniProt Accession	P24807
Predicted Molecular Weight	8kDa
Shipping Condition	Ice bag