

Rabbit Anti-Mouse Cd24 Monoclonal Antibody For If/Icc And Fc - P24807

Item Number: CM0005690

Introduction

-conjugated rabbit monoclonal antibody against mouse CD24 for immunofluorescence (IF/ICC) and flow cytometry (FC). CD24 is a GPI-anchored cell surface glycoprotein involved in cell differentiation, thymocyte development, B-cell activation, and autoimmunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
GeneID (Human)	100133941
GeneID	12484
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag