

## Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Fc - P01590

Item Number: CM0008519

### Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD25 (IL-2 receptor alpha). Specifically reacts with mouse CD25. Validated for flow cytometry applications in immunology and regulatory T cell research, suitable for studying immune tolerance mechanisms.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD25 mAb                                                                                                            |
| Remarks/Alias    | CD25; IL2r; Ly-43                                                                                                                     |
| Species          | Mouse                                                                                                                                 |
| Gene ID (Human)  | 3559                                                                                                                                  |
| Gene ID          | 16184                                                                                                                                 |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | FC                                                                                                                                    |
| Cross Reactivity | Mouse                                                                                                                                 |
| SWISS            | P01590                                                                                                                                |
| Protein Weight   | 30KDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |