

Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Flow Cytometry - P01590

Item Number: CM0008908

Introduction

-conjugated rabbit monoclonal antibody against mouse CD25 (IL-2R α) for flow cytometry. Validated in mouse. Targets interleukin-2 receptor subunit alpha, a key regulator of Treg-mediated immune tolerance.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
GeneID (Human)	3559
GeneID	16184
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25. (NP_032393.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	30KDa
Shipping	Ice bag