

Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Flow Cytometry - P01590

Item Number: CM0011348



Introduction

High-quality ABflo 450-conjugated rabbit anti-mouse CD25 monoclonal antibody for flow cytometry. Recognizes mouse CD25 (IL-2 receptor alpha), enabling Treg identification and immune tolerance research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
GenelD (Human)	3559
GenelD	16184
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	31kDa
Shipping	Ice bag