

Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Flow Cytometry - P01590

Item Number: CM0007847

Introduction

Rabbit monoclonal antibody targeting mouse CD25 (IL2RA), validated for flow cytometry. Useful for Treg identification and IL-2 receptor studies. Specific to mouse, 30kDa. - P01590

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
GeneID (Human)	3559
GeneID	16184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	30kDa
Shipping	Ice bag