

## Rabbit Anti-Mouse Cd25 Monoclonal Antibody (Il2Ra, P55) For Flow Cytometry - P01590

Item Number: CM0012578

### Introduction

Rabbit anti-Mouse CD25 monoclonal antibody, validated for flow cytometry. Targets mouse IL2RA/p55 (UniProt P01590, 30 kDa). Ideal for Treg and immune tolerance flow cytometry research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD25 mAb                                                                                       |
| Remarks/Alias    | CD25; Il2r; Ly-43                                                                                                |
| Species          | Mouse                                                                                                            |
| Gene ID (Human)  | 3559                                                                                                             |
| Gene ID          | 16184                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3) |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | FC                                                                                                               |
| Cross Reactivity | Mouse                                                                                                            |
| SWISS            | P01590                                                                                                           |
| Protein Weight   | 30kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |