


Rabbit Anti-Mouse Cd252 (Ox40 Ligand) Monoclonal Antibody For If/Icc, Fc - P43488

Item Number: CM0009979



Introduction

-conjugated rabbit monoclonal antibody specific for mouse CD252 (OX40 Ligand). Validated for IF/ICC and FC. Co-stimulates T-cell proliferation and cytokine production. ~22 kDa membrane protein.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD252 (OX40 Ligand) mAb
Remarks/Alias	Ath1; gp34; Ath-1; Ox40L; TXGP1; CD134L; OX-40L; Tnlg2b; Txgp1l
Species	Mouse
Gene ID (Human)	7292
Gene ID	22164
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 51-198 of mouse CD252 (OX40 Ligand) (NP_033478.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P43488
Protein Weight	22kDa
Shipping	Ice bag