

Rabbit Anti-Mouse Cd262/Trailr2 Monoclonal Antibody For Fc - Q9Qzm4

Item Number: CM0010289

Introduction

High-quality -conjugated rabbit monoclonal antibody specific for mouse CD262/TRAILR2 (Death receptor 5). Validated for flow cytometry to investigate apoptosis, TRAIL signaling, and ER stress in heart, lung, and kidney tissues.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD262/TRAILR2 mAb
Remarks/Alias	MK; DR5; Ly98; KILLER; TRICKB; TRAILR2; TRICK2A; TRICK2B
Species	Mouse
Gene ID (Human)	8795
Gene ID	21933
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 53-170 of mouse CD262/TRAILR2(NP_064671.2).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QZM4
Protein Weight	42kDa
Shipping	Ice bag