

Rabbit Anti-Mouse Cd262/Trailr2 Monoclonal Antibody For Fc - Q9Qzm4

Item Number: CM0008404

Introduction

-conjugated rabbit monoclonal antibody against mouse CD262/TRAILR2 (Death receptor 5, MK). Validated for flow cytometry. Detects a 42 kDa protein, reacts with mouse. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD262/TRAILR2 mAb
Remarks/Alias	MK; DR5; Ly98; KILLER; TRICKB; TRAILR2; TRICK2A; TRICK2B
Species	Mouse
GenelD (human)	8795
GenelD	21933
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 53-170 of mouse CD262/TRAILR2 (NP_064671.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QZM4
Protein Weight	42kDa
Shipping	Ice bag