

Rabbit Anti-Mouse Cd262/Trailr2 Monoclonal Antibody For Flow Cytometry - Q9Qzm4

Item Number: CM0005340

Introduction

-conjugated rabbit monoclonal antibody against mouse CD262/TRAILR2 (DR5), a death receptor that mediates apoptosis upon TRAIL binding. Designed for flow cytometry detection in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD262/TRAILR2 mAb
Remarks/Alias	MK, DR5, Ly98, KILLER, TRICKB, TRAILR2, TRICK2A, TRICK2B
Species	Mouse
GeneID (Human)	8795
GeneID	21933
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 53-170 of mouse CD262/TRAILR2 (NP_064671.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QZM4
Protein Weight	42kDa
Shipping	Ice bag