

Rabbit Anti-Mouse Cd267/Taci Monoclonal Antibody For Fc - Q9Et35

Item Number: CM0011377

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD267/TACI (TNFRSF13B). Validated in flow cytometry (FC). Reacts with mouse. Useful for B-cell and T-cell function studies, humoral immunity research, and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD267/TACI mAb
Remarks/Alias	Taci; 1200009E08Rik
Species	Mouse
Gene ID (Human)	23495
Gene ID (Mouse)	57916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 5-129 of mouse CD267/TACI (NP_067324.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9ET35
Protein Weight	27kDa
Shipping	Ice bag