

Abflo 647 Rabbit Anti-Mouse Cd267/Taci Monoclonal Antibody For Flow Cytometry - Q9Et35

Item Number: CM0012452

Introduction

ABflo 647-conjugated rabbit monoclonal antibody against mouse CD267/TACI, validated for flow cytometry. Targets a receptor involved in B and T cell immunity and humoral regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD267/TACI mAb
Remarks/Alias	Taci; 1200009E08Rik
Species	Mouse
Gene ID (Human)	23495
Gene ID	57916
Immunogen	Recombinant Protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9ET35
Protein Weight	27kDa
Shipping	Ice bag