

Rabbit Anti-Mouse Cd270/Hvem Monoclonal Antibody For Fc - Q80Wm9

Item Number: CM0008893



Introduction

-conjugated rabbit monoclonal antibody specific for mouse CD270/HVEM. Validated for flow cytometry (FC), enabling detection of this immune checkpoint receptor in research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD270/HVEM mAb
Remarks/Alias	TR2; Atar; HveA; Hvem; Tnfrs14
Species	Mouse
GeneID (Human)	8764
GeneID	230979
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 39-207 of mouse CD270/HVEM (NP_849262.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q80WM9
Protein Weight	30kDa
Shipping	Ice bag