

Rabbit Anti-Mouse Cd272/Btla Mab For Flow Cytometry - Q7Tsa3

Item Number: CM0005888



Introduction

-conjugated rabbit monoclonal antibody specific for mouse CD272/BTLA, validated for flow cytometry. Suitable for studying inhibitory receptor function and immune regulation in T and B cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD272/BTLA mAb
Remarks/Alias	A630002H24
Species	Mouse
Gene ID (Human)	151888
Gene ID	208154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-176 of mouse CD272/BTLA (NP_001032808.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q7TSA3
Protein Weight	34kDa
Shipping	Ice bag