


Rabbit Anti-Mouse Cd274/Pd-L1 Mab For Flow Cytometry - Q9Ep73

Item Number: CM0012507



Introduction

-conjugated rabbit anti-mouse CD274/PD-L1 monoclonal antibody validated for flow cytometry. Detects PD-L1 (B7-H1, Q9EP73) on mouse cells. Ideal for immuno-oncology studies. Supplied by CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD274/PD-L1 mAb
Remarks/Alias	B7h1; Pdl1; Pdcd1l1; Pdcd1lg1; A530045L16Rik
Species	Mouse
Gene ID (Human)	29126
Gene ID	60533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-235 of mouse CD274/PD-L1(NP_068693.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9EP73
Protein Weight	32kDa
Shipping	Ice bag