


Rabbit Anti-Mouse Cd274/Pd-L1 Monoclonal Antibody For Flow Cytometry - Q9Ep73

Item Number: CM0010193



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD274/PD-L1 for flow cytometry applications. Reacts with mouse samples; 32 kDa. Ideal for studying immune checkpoint pathways and anti-tumor immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD274/PD-L1 mAb
Remarks/Alias	B7h1, Pdl1, Pdcd1l1, Pdcd1lg1, A530045L16Rik
Species	Mouse
GeneID (Human)	29126
GeneID	60533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-235 of mouse CD274/PD-L1 (NP_068693.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9EP73
Protein Weight	32kDa
Shipping	Ice bag