

Rabbit Anti-Mouse Cd28 Monoclonal Antibody For Flow Cytometry - P31041

Item Number: CM0012097



Introduction

-conjugated rabbit monoclonal antibody against mouse CD28 for flow cytometry. Specific for mouse CD28, recognizes a 25 kDa protein. Suitable for T-cell immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD28 mAb
Remarks/Alias	Cd28
Species	Mouse
Gene ID (Human)	940
Gene ID	12487
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-149 of mouse CD28 (NP_031668.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P31041
Protein Weight	25kDa
Shipping	Ice bag