

Rabbit Anti-Mouse Cd301 (Mgl1/2) Monoclonal Antibody For Fc - P49300 / Q8Jzn1

Item Number: CM0010921

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD301 (MGL1/2), suitable for flow cytometry. Recognizes Clec10a/Clec10b, key macrophage lectins in immunity and metabolism.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD301 (MGL1/2) mAb
Remarks/Alias	Mgl; Mgl1; CD301a; M-ASGP-BP-1
Species	Mouse
Gene ID (Human)	10462/54551
Gene ID	17312/216864
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 57-304 of mouse CD301a (NP_034926.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P49300/Q8JZN1
Protein Weight	34kDa/38kDa
Shipping	Ice bag