

Rabbit Anti-Mouse Cd301 (Mgl1/2) Monoclonal Antibody For Flow Cytometry - P49300 / Q8Jzn1

Item Number: CM0011628

Introduction

-conjugated rabbit monoclonal antibody against mouse CD301 (MGL1/2). Specific for flow cytometry, reacts with mouse samples. Recognizes macrophage galactose-type lectins, useful for immunology and infection studies.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD301 (MGL1/2) mAb
Remarks/Alias	Mgl; Mgl1; CD301a; M-ASGP-BP-1
Species	Mouse
GeneID (Human)	10462/54551
GeneID	17312/216864
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 57-304 of mouse CD301a (NP_034926.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P49300/Q8JZN1
Protein Weight	34kDa/38kDa
Shipping	Ice bag