

Rabbit Anti-Mouse Cd301 (Mgl1/2) Monoclonal Antibody For Flow Cytometry - P49300 / Q8Jzn1

Item Number: CM0010078

Introduction

-conjugated rabbit anti-mouse CD301 (MGL1/2) monoclonal antibody for flow cytometry. Detects mouse CD301a and CD301b, galactose/N-acetylgalactosamine-specific lectins on activated macrophages, involved in immune responses. Validated in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD301 (MGL1/2) mAb
Remarks/Alias	Mgl; Mgl1; CD301a; M-ASGP-BP-1
Species	Mouse
GeneID (Human)	10462/54551
GeneID	17312/216864
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 57-304 of mouse CD301a (NP_034926.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P49300/Q8JZN1
Protein Weight	34kDa/38kDa
Shipping	Ice bag