


Rabbit Anti-Mouse Cd317/Bst2 Monoclonal Antibody For Flow Cytometry - Q8R2Q8

Item Number: CM0010536



Introduction

Rabbit monoclonal antibody targeting mouse CD317/BST2 (GREG, HM1.24), validated for flow cytometry. Recognizes a 19kDa bone marrow stromal antigen involved in antiviral restriction. Suitable for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD317/BST2 mAb
Remarks/Alias	GREG; Bst-2; CD317; DAMP-1; 2310015I10Rik
Species	Mouse
GeneID (Human)	684
GeneID	69550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-152 of mouse CD317/BST2 (NP_932763.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8R2Q8
Protein Weight	19kDa
Shipping	Ice bag