

Rabbit Anti-Mouse Cd34 Mab - Q64314

Item Number: CM0010351

Introduction

Rabbit anti-Mouse CD34 monoclonal antibody for IF/ICC and FC. Targets hematopoietic progenitor cell antigen CD34 (Q64314). Ideal for mouse stem cell research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Mouse
GeneID (human)	947
GeneID	12490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-287 of mouse CD34 (NP_598415.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	35kDa/41kDa
Shipping	Ice bag