

Rabbit Anti-Mouse Cd34 Monoclonal Antibody For Fc - Q64314

Item Number: CM0004836

Introduction

Rabbit anti-Mouse CD34 monoclonal antibody for flow cytometry. Recognizes mouse CD34 (41kDa), involved in hematopoietic progenitor cell adhesion and stem cell attachment. Suitable for mouse FC assays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD34 mAb
Remarks/Alias	CD34;CD34 molecule;GIG3;MORT1
Species	Mouse
GeneID (Human)	947
GeneID	12490
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 35-287 of mouse CD34 (NP_598415.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	41kDa
Shipping	Ice bag