

Rabbit Anti-Mouse Cd34 Monoclonal Antibody For Flow Cytometry - Q64314

Item Number: CM0009441

Introduction

Rabbit monoclonal antibody against mouse CD34, conjugated with for flow cytometry. Targets hematopoietic progenitor cell antigen CD34, a key marker in stem cell and hematopoiesis research. Suitable for mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Mouse
Gene ID (Human)	947
Gene ID	12490
Immunogen	Recombinant protein (aa 35-287 of mouse CD34, NP_598415.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	41 kDa
Shipping	Ice bag