

Rabbit Anti-Mouse Cd352/Slamf6 Monoclonal Antibody For Flow Cytometry - Q9Et39

Item Number: CM0009119

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD352/SLAMF6 for flow cytometry. Suitable for studying immune cell modulation, NK cell activation, and T-cell differentiation in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD352/SLAMF6 mAb
Remarks/Alias	KAL1; NTBA; KAL1b; Ly108; NTB-A; SF2000
Species	Mouse
GenelD (Human)	114836
GenelD (Mouse)	30925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-235 of mouse CD352/SLAMF6 (NP_001334115.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9ET39
Protein Weight	39kDa
Shipping	Ice bag