

Rabbit Anti-Mouse Cd366 (Tim-3) Mab For Flow Cytometry - Q8Vim0

Item Number: CM0010198



Introduction

-conjugated rabbit monoclonal antibody against mouse TIM-3 (CD366) for flow cytometry. Specific for mouse, no cross-reactivity. Suitable for studying immune checkpoint regulation in Th1, CD8+ T cells, NK cells, macrophages.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD366 (TIM-3) mAb
Remarks/Alias	Tim3; TIM-3; Timd3
Species	Mouse
GeneID (Human)	84868
GeneID	171285
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8VIM0
Protein Weight	22kDa/31kDa
Shipping	Ice bag