

Rabbit Anti-Mouse Cd38 Monoclonal Antibody For Flow Cytometry - P56528

Item Number: CM0007922



Introduction

-conjugated rabbit anti-mouse CD38 monoclonal antibody for flow cytometry. Detects CD38, a multifunctional transmembrane glycoprotein involved in calcium mobilization, lymphocyte activation, and receptor-mediated signaling. Suitable for immune cell profiling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD38 mAb
Remarks/Alias	I-19; ADPRC 1; Cd38-rs1
Species	Mouse
GeneID (human)	952
GeneID	12494
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 45-304 of mouse CD38 (NP_031672.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P56528
Protein Weight	34kDa
Shipping	Ice bag