

Rabbit Anti-Mouse Cd4 Mab For Fc - P06332

Item Number: CM0005847

Introduction

Rabbit anti-Mouse CD4 monoclonal antibody validated for flow cytometry. Recognizes mouse CD4 (L3T4), a T-helper cell coreceptor. Suitable for immunophenotyping and immune response research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
GeneID (Human)	920
GeneID	12504
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag