

Rabbit Anti-Mouse Cd4 Monoclonal Antibody For Flow Cytometry - P06332

Item Number: CM0009740

Introduction

-conjugated rabbit anti-mouse CD4 monoclonal antibody for flow cytometry. Validated for mouse samples. Ideal for T-helper cell identification and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
GeneID (Human)	920
GeneID (Mouse)	12504
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag