

Abflo 488 Rabbit Anti-Mouse Cd4 Monoclonal Antibody For Flow Cytometry - P06332

Item Number: CM0005612

Introduction

Rabbit monoclonal antibody against mouse CD4, conjugated with ABflo 488, specifically designed for flow cytometry. Recognizes mouse CD4 (51 kDa), a T-helper cell marker. Suitable for immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
Gene ID (Human)	920
Gene ID (Mouse)	12504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag