

Rabbit Anti-Mouse Cd40 Monoclonal Antibody For Flow Cytometry - P27512

Item Number: CM0011792

Introduction

-conjugated rabbit anti-mouse CD40 monoclonal antibody validated for flow cytometry. Targets mouse CD40 (UniProt P27512), a TNF receptor superfamily member involved in B cell and macrophage activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD40 mAb
Remarks/Alias	IGM; p50; Bp50; GP39; IMD3; TRAP; HIGM1; T-BAM; Tnfrsf5
Species	Mouse
Gene ID (Human)	958
Gene ID	21939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-189 of mouse CD40 (NP_035741.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P27512
Protein Weight	32kDa
Shipping	Ice bag