


Rabbit Anti-Mouse Cd44 Monoclonal Antibody For Flow Cytometry - P15379

Item Number: CM0012292



Introduction

-conjugated rabbit monoclonal antibody against mouse CD44 (Ly-24/Pgp-1). Designed for flow cytometry, recognizes native mouse CD44 antigen. CD44 plays critical roles in cell adhesion, lymphocyte homing, and tumor metastasis. Ideal for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
GeneID (Human)	960
GeneID	12505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	85kDa
Shipping	Ice bag