


## Rabbit Anti-Mouse Cd44 Monoclonal Antibody ( ) For Flow Cytometry - P15379

Item Number: CM0011806



### Introduction

-conjugated rabbit monoclonal antibody against mouse CD44 antigen, validated for flow cytometry. Targets hyaluronate receptor involved in lymphocyte homing, cell adhesion, and migration.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD44 mAb                                                                                        |
| Remarks/Alias    | Ly-24; Pgp-1; HERMES                                                                                              |
| Species          | Mouse                                                                                                             |
| Gene ID (Human)  | 960                                                                                                               |
| Gene ID (Mouse)  | 12505                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2). |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Mouse                                                                                                             |
| SWISS            | P15379                                                                                                            |
| Protein Weight   | 85kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |