

Rabbit Anti-Mouse Cd44 Monoclonal Antibody For If/Icc, Fc - P15379

Item Number: CM0008848



Introduction

-conjugated rabbit monoclonal antibody against mouse CD44 (Ly-24/Pgp-1/HERMES), validated for immunofluorescence (IF/ICC) and flow cytometry (FC). Suitable for studying cell adhesion, migration, and lymphocyte homing in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD44 mAb
Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
Gene ID (Human)	960
Gene ID	12505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	85kDa
Shipping	Ice bag