

## Rabbit Anti-Mouse Cd45 Mab - P06800

Item Number: CM0008979

### Introduction

Rabbit anti-Mouse CD45 monoclonal antibody for flow cytometry. Specific to mouse CD45, validated for FC, with protein weights of 130/140/144 kDa. Ideal for immune cell identification and research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD45 mAb                                                                                          |
| Remarks/Alias    | loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyl-4                                                                     |
| Species          | Mouse                                                                                                               |
| GeneID (Human)   | 5788                                                                                                                |
| GeneID (Mouse)   | 19264                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | FC                                                                                                                  |
| Cross Reactivity | Mouse                                                                                                               |
| SWISS            | P06800                                                                                                              |
| Protein Weight   | 130kDa/140kDa/144kDa                                                                                                |
| Shipping         | Ice bag                                                                                                             |