

Rabbit Anti-Mouse Cd45 Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0011125

Introduction

-conjugated rabbit monoclonal antibody against mouse CD45, validated for flow cytometry. Recognizes isoforms ranging from 130 to 144 kDa. Ideal for immune cell profiling and tyrosine phosphatase signaling research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
Gene ID (Human)	5788
Gene ID	19264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag